

# X-RAY GENERATORS

# CMP 200<sup>®</sup>

(SERIES)



**CPI**

## CMP 200<sup>®</sup> series

True 100 kHz + (RAD only)  
'clinical-quality' X-Ray generators  
that promise the best price/  
performance ratio in their market.

### 3 different models to choose from:

- 32 kW and 40 kW
  - 1 Phase 208/230 VAC (125kV) only
  - 3 Phase 208/230 VAC (125 kV) only
  - 3 Phase 400/480 (125 kV) standard  
(150 kV optional)
- 50 kW
  - 3 Phase 400/480 VAC (150 kV) standard

Cost-effective for all radiographic  
applications (hospitals, medical clinics,  
chiropractic and veterinary applications)

**Small. Compact. Lightweight.**  
A low cost alternative.

X-RAY  
**Generators  
of Choice**



# Family of X-Ray 'Generators-of-Choice'.



## Look at these Features

- State-of-the-art 100kHz + FET-based switching technology from the Indico 100® is also here in the CMP 200®.
- Anatomical programs available in several languages – and easily edited using the Operator Console
- Technique Selection:
  1. kV with AEC (Optional)
  2. kV/mAs
  3. kV/mA/Time
  4. kV/cm Thickness
- Low kV ripple -- maximizes image quality by providing constant potential performance at a fraction of the cost of more expensive models.
- Patient Body-Part Thickness mode
- CPI's output frequencies in excess of 200 kHz are the highest in the industry and produce the lowest kV ripple, resulting in the highest potential image quality and the lowest stored energy. This simply means you get maximum tube protection and the longest tube life.
- Operator Control Console available with international symbols or text.
- Touchscreen Control Console (Optional)
- Wall Mount (Optional)
- Automatic Exposure Control (AEC) (Optional)
- Dose Area Product (DAP) Interface (Optional)
- Automatic Tube Calibration
- Compact, cost effective design saves you money in transportation and installation.
- Handswitch (Optional)
- Collimator power supply
- Tomography interface (Optional)



## World Class Product Support

Our greatest strength is our world class reputation for customer support. A dedicated support team is only a phone call away at 905-877-0161 – or at our North American toll-free number 1-888-CPI-XRAY (274-9729)

Call us. We are here to help you with all of your installation needs.



Communications & Power Industries Canada Inc.  
communications & medical products division  
45 River Drive / Georgetown, Ontario / Canada L7G 2J4  
Hot Line Telephone: 1-888-CPI-XRAY  
TEL: 905-877-0161 / FAX: 905-877-5327  
E-MAIL: [marketing@cmp.cpii.com](mailto:marketing@cmp.cpii.com)  
WEB: [www.cpii.com/cmp](http://www.cpii.com/cmp)



Some features shown may be optional. Please refer to Data Sheet.